


WELLMEDPRO

Hospital Furnitures



Contents :

001	-----	Electric 5-Function Patient Bed
002	-----	Electric 3-Function & Manual 3-Function Hospital Bed
003	-----	Manual 3-Crank / 2-Crank / Single-Crank Patient Bed
004	-----	Patient Transfer Cart
005	-----	Stretcher & Wheel Chair
006	-----	Examination Couch / Bedside Cabinet / Overbed Dining Table / Chairs / Wards Screen / Drip Stand
007	-----	Waiting Chair / Infusion Chair / Instrument Trolley / Surgeon Stool / Kidney Tray / Sterilization Drum / Stainless Sink
008	-----	Medical Trolley / Medical Waste Bin

WELLMEDPRO

Wellview Universal Ltd.

Address: Room 1617-1618 Star House, 3 Salisbury Road,
TST, Kowloon, Hong Kong

E-mail: nacho@wellviewuniversal.net
www.wellmedpro.com



Hospital Furnitures

Electric Patient Bed



YSHB-LZ5D
Double Column Motors

- ★ Main functions: A: Back section up/down;
B: Knee section up/down;
C: Whole bed up/down;
D: Trendelenburg and Reverse Trendelenburg;
- ★ Overall dimension: L2140xW1050xH505-805mm
- ★ Back rest tilting: 75°±5°; Knee board tilting: 35°±5°;
- ★ Hi-lo adjustment between: 505-805mm;
- ★ Reverse and trendelenburg 15°±2°;
- ★ Standard configuration: 2 pcs column motors; 2 pcs advanced medical actuators; 1 pc control box; weighing function; X-ray transmission; 1 pc ACP(Attendant Control Operation Panel); Central-controlled brake system; 4pcs luxurious centrallocking castors; 1 pair plastic head/foot-end board; 1 pair Luxurious PP siderails; 4 pcs urine hooks; 1pc S.S I.V pole; CPR function; 1 pc emergency stop switch; 1 pc Backup power system (UPS);



YSHB-D504
Weighing Function

- ★ Main functions: A: Back section up/down.
B: Knee section up/down.
C: Height adjustment.
D: Trendelenburg and Reverse Trendelenburg.
- ★ Overall dimension: L2220xW1050xH490-795mm±10mm
- ★ Back rest tilting: 75°±5°; thigh rest tilting: 40°±5°;
- ★ Hi-lo adjustment between: 490-795mm±10mm;
- ★ Reverse and trendelenburg:12°±2°;
- ★ Mattress platform:1930x890mm;
- ★ Weight:153KG; Package size: 2070mmx980mmx490mm
- ★ Standard configuration: 4pcs medical actuators; 1pc control box; 1pc nurse console; weighing system function; central-controlled brake system, 4pcs luxurious central locking castors; 1 pair plastic head/foot-end board; 1 pair PP siderails; 4pcs urine hooks; 1pc S.S I.V pole; 1pc Backup power system(UPS); CPR function; underbed LED light; 2pcs angle indicator: 0-90°.



YSHB-D503
Standing Bed

- ★ Main functions: A: Back section up/down;
B: Knee section up/down;
C: Whole bed up/down;
D: Trendelenburg and Reverse Trendelenburg;
E: Standing function;
- ★ Overall dimension: L2250xW1080xH490-880mm
- ★ Back rest tilting: 70°±5°; Knee board tilting: 30°±5°;
- ★ Hi-lo adjustment between: 490-880mm;
- ★ Reverse and trendelenburg: 15°±2°;
- ★ Standing functional angle: 85°±5°;
- ★ Standard configuration:
5 pcs Advanced medical actuators,
1 pc control box, 1pc handset;
Central-controlled brake system,
6 pcs luxurious central locking castors; 1 pair plastic head/foot-end boards;
1 pair Luxurious PP siderails;
4 pcs urine hooks; 1 pc S.S I.V pole;
- ★ Options: Backup power system (UPS); CPR function; Medical mattress; Dining board; Bedside cabinet;



YSHB-D501

- ★ Main functions: A: Back section up/down.
B: Knee section up/down.
C: Height adjustment.
D: Trendelenburg and Reverse Trendelenburg.
- ★ Overall dimension: L2170xW990xH460-720mm±10mm
- ★ Back rest tilting:75°±5°; thigh rest tilting: 40°±5°;
- ★ Hi-lo adjustment between:460-720mm±10mm;
- ★ Reverse and trendelenburg: 15°±2°.
- ★ Standard configuration: 4pcs medical actuators; 1pc control box; 1pc handset; central-controlled brake system, 4pcs luxurious central locking castors; 1 pair plastic head/foot-end board; 1 pair PP siderails; 4pcs urine hooks; 1pc S.S I.V pole;
- ★ Options: backup power system(UPS); CPR function; medical mattress; overbed table; bedside cabinet; patient file panel.



YSHB-HN05E

- ★ Main functions: A: Back section up/down.
B: Knee section up/down.
C: Height adjustment.
D: Trendelenburg and Reverse Trendelenburg.
- ★ Overall dimension: L2050xW1050xH420-680mm±10mm
- ★ Back rest tilting:75°±5°; thigh rest tilting: 35°±5°;
- ★ Hi-lo adjustment between:420-680mm±10mm;
- ★ Reverse and trendelenburg: 12°±2°.
- ★ Standard configuration: motor; 1pc control box; 1pc handset; brake system, 4pcs castors; 1 pair plastic head/foot-end board; 1 pair PP siderails; 1pc S.S I.V pole;
- ★ Options: medical mattress; overbed table; bedside cabinet; patient file panel.



YSHB-HN05D

- ★ Main functions: A: Back section up/down.
B: Knee section up/down.
C: Height adjustment.
D: Trendelenburg and Reverse Trendelenburg.
- ★ Overall dimension: L2050xW960xH420-680mm±10mm
- ★ Back rest tilting:75°±5°; thigh rest tilting: 35°±5°;
- ★ Hi-lo adjustment between:420-680mm±10mm;
- ★ Reverse and trendelenburg: 12°±2°.
- ★ Standard configuration: motor; 1pc control box; 1pc handset; brake system, 4pcs castors; 1 pair plastic head/foot-end board; 1 pair aluminum alloy siderails; 1pc S.S I.V pole;
- ★ Options: medical mattress; overbed table; bedside cabinet; patient file panel.



Hospital Furnitures

Electric & Manual Patient Bed



YSHB-D303
Electric 3-Function

- ★ Main functions: A: Back section up/down.
B: Knee section up/down.
C: Height adjustment
- ★ Overall dimension : L2170*W990*H460-720mm±10mm
- ★ Back rest tilting :80°±5°, thigh rest tilting: 40°±5°
- ★ Hi-lo adjustment between: 460-720mm±10mm;
- ★ Standard configuration:
3pcs medical actuators, 1pcs control box, 1pcs handset, central-controlled brake system, 4pcs central locking castors; 1 pair ABS head/foot-end board; 1 pair PP siderails; 4pcs urine hooks; 1pcs S.S I.V pole;
- ★ Options: Backup power system (UPS); CPR function; medical mattress, dining board; bedside cabinet; patient file panel ; sundries basket



YSHB-HN03E
Electric 3-Function

- ★ Main functions: A: Back section up/down.
B: Knee section up/down.
C: Height adjustment.
- ★ Overall dimension: L2115xW960xH420-680mm±10mm
- ★ Back rest tilting: 75°±5°; thigh rest tilting: 35°±5°;
- ★ Hi-lo adjustment between: 420-680mm±10mm;
- ★ Standard configuration:
3pcs electric motor;
1pc control box;
1pc handset;
brake system,
4pcs castors;
1 pair ABS head/foot-end board;
1 pair PP siderails;
1pc S.S I.V pole;
- ★ Options: medical mattress; overbed table; bedside cabinet; patient file panel.



YSHB-HN03D
Electric 3-Function

- ★ Main functions: A: Back section up/down.
B: Knee section up/down.
C: Height adjustment.
- ★ Overall dimension: L2115xW960xH420-680mm±10mm
- ★ Back rest tilting: 75°±5°; thigh rest tilting: 35°±5°;
- ★ Hi-lo adjustment between: 420-680mm±10mm;
- ★ Standard configuration:
3pcs electric motors;
1pc control box;
1pc handset; brake system,
4pcs castors;
1 pair ABS head/foot-end board;
1 pair aluminum alloy siderails;
1pc S.S I.V pole;
- ★ Options: medical mattress; overbed table; bedside cabinet; patient file panel.



YSHB-HN03B
Manual 3-Crank

- ★ Main functions: A: Back section up/down.
B: Knee section up/down.
C: Height adjustment.
- ★ Overall dimension: L2115xW960xH420-680mm±10mm
- ★ Back rest tilting: 75°±5°; thigh rest tilting: 35°±5°;
- ★ Hi-lo adjustment between: 420-680mm±10mm;
- ★ Standard configuration:
4pcs castors;
1 pair ABS head/foot-end board;
1 pair PP siderails;
1pc S.S I.V pole;
- ★ Options: medical mattress; overbed table; bedside cabinet; patient file panel.



YSHB-HN03A
Manual 3-Crank

- ★ Main functions: A: Back section up/down.
B: Knee section up/down.
C: Height adjustment.
- ★ Overall dimension: L2115xW960xH420-680mm±10mm
- ★ Back rest tilting: 75°±5°; thigh rest tilting: 35°±5°;
- ★ Hi-lo adjustment between: 420-680mm±10mm;
- ★ Standard configuration:
4pcs castors;
1 pair ABS head/foot-end board;
1 pair aluminum alloy siderails;
1pc S.S I.V pole;
- ★ Options: medical mattress; overbed table; bedside cabinet; patient file panel.



YSGH1015-B
Orthopedic Traction Bed

- ★ Main functions: A: Back section up/down.
B: Left foot up/down.
C: Right foot up/down.
D: Height adjustment.
- ★ Overall dimension: L2115xW980xH500-2000mm±10mm
- ★ Back rest tilting: 80°±5°; thigh rest tilting: 40°±5°;
- ★ Standard configuration:
1 pair ABS head/foot-end board;
1 pair aluminum alloy siderails;
1pc S.S I.V pole;
1 set orthopedic rack
- ★ Options: medical mattress, central control castors, overbed table, bedside cabinet



Hospital Furnitures

Manual Patient Bed



YSHB-S231



- ★ Main functions: A: Back section up/down.
B: Knee section up/down.
- ★ Overall dimension: L2100*W970*H480mm±10mm
- ★ Back rest tilting: 80°±5°, thigh rest tilting: 40°±5°
- ★ Standard configuration:
 - 2 sets manual crank mechanism systems;
 - 4pcs castors;
 - 1 pair ABS head/foot-end boards;
 - 1 pair Aluminium siderails;
 - 4pcs urine hooks;
 - 1pc S.S I.V pole;
- ★ Options: Medical mattress, dining board, bedside cabinet, sundries basket.



YSHB-HN02B



- ★ Main functions: A: Back section up/down.
B: Knee section up/down.
- ★ Overall dimension: L2115xW960xH500mm±10mm
- ★ Back rest tilting: 75°±5°; thigh rest tilting: 35°±5°;
- ★ Standard configuration:
 - 4pcs castors;
 - 1 pair ABS head/foot-end board;
 - 1 pair aluminum alloy siderails;
 - 1pc S.S I.V pole;
- ★ Options: medical mattress; overbed table; bedside cabinet; patient file panel.



YSHB-HN02A



- ★ Main functions: A: Back section up/down.
B: Knee section up/down.
- ★ Overall dimension: L2115xW960xH500mm±10mm
- ★ Back rest tilting: 75°±5°; thigh rest tilting: 35°±5°;
- ★ Standard configuration:
 - 4pcs castors;
 - 1 pair ABS head/foot-end board;
 - 1 pair aluminum alloy siderails;
 - 1pc S.S I.V pole;
- ★ Options: medical mattress; overbed table; bedside cabinet; patient file panel.



YSHB-HN01B



- ★ Main functions: A: Back section up/down.
- ★ Overall dimension: L2115xW960xH500mm±10mm
- ★ Back rest tilting: 75°±5°
- ★ Standard configuration:
 - 4pcs castors;
 - 1 pair ABS head/foot-end board;
 - 1 pair aluminum alloy siderails;
 - 1pc S.S I.V pole;
- ★ Options: medical mattress; overbed table; bedside cabinet; patient file panel.



YSHB-HN01A



- ★ Main functions: A: Back section up/down.
- ★ Overall dimension: L2115xW960xH500mm±10mm
- ★ Back rest tilting: 75°±5°
- ★ Standard configuration:
 - 4pcs castors;
 - 1 pair ABS head/foot-end board;
 - 1 pair aluminum alloy siderails;
 - 1pc S.S I.V pole;
- ★ Options: medical mattress; overbed table; bedside cabinet; patient file panel.



YSGH4901 Pediatric Bed



- ★ Size (mm): 1600*700*550
- ★ Materials: Coated steel
- ★ Colors: White
- ★ Standard: Side rails/IV Pole
- ★ Optionals: Castors/Mattress



Hospital Furnitures

Patient Transfer Cart



YSHB-KX857



- ★ Hand rescue vehicles, hand shake clockwise, counter clock wise direction to the body, shaking, body down;
- ★ A clear arrow mark indicating the user to operate
- ★ In the control wheel, diameter 150mm, equipped with a central steering wheel to facilitate the user to control the direction of the vehicle body;
- ★ High-low range: 820mm - 520mm
- ★ Standare accessories: mattress, I.V.Pole frame;
- ★ Dimensions: 1930×650×(820mm-520mm)



YSHB-KX882



- ★ Dimension: 2130×750mm
- ★ Two Separate hydraulic pump, bed position adjustable.
- ★ Back section panel and Leg section panel can be adjustable, Trendelenburg position, Anti-Trendelenburg position, Back-rest can through X-ray.
- ★ Central lock brake, diameter 150mm, 6 inch, mate with guide wheel, operate the direction of cart.
- ★ Top Coat: Electrostatic powder coating; Bed Frame: Cold-rolled steel tube.
- ★ Optional: Pemissible X-ray Stainless steel and aluminum alloy side rail.
- ★ Accessories: I.V.Pole, Mattress



YSHB-KX896



- ★ It is a manual connecting transfer stretcher, which can be used to avoid cross-infection
- ★ Cart A inside the operation room, Cart B outside, adjusting them to be the same height, then push Cart B to Cart A, two carts will be connected, following the orbit, bed board C can be pushed from one cart to another, and VICE VERSA;
- ★ Bright arrow symbols direct the users how to operate clearly;
- ★ Hi-low range: 920mm-645mm,
- ★ Connecting transfer Stretcher, exterior dimension: 3850×800mm, bank angle display: 0~60°±5°;
- ★ 4 pieces of lifting PP guardrail, integral injection molding;
- ★ Standard accessories: mattress, I.V.pole.



YSHB-ES005



- ★ Main functions: A.Back rest up/down
B.Thigh rest up/down
C.Height adjustment
D.Trendelenburg and revise trendelenburg
- ★ Size: L2080*W850*H700-1000mm±10mm
- ★ Back rest tilting : 70°±5°
- ★ Thigh rest tilting : 30°±5°
- ★ Trendelenburg : 15°±2°
- ★ Reverse trendelenburg:15°±2°
- ★ Standard configuration:
 - 1 set manual crank mechanism system
 - 2 back rise/ fall gas springs
 - 2 rise and fall hydraulic colomnns
 - 1 pair steel side rails(6 bars each piece)
 - 4 pcs central locking castors of 6 inch
 - 1 directional wheel (5th wheel)
 - 1 oxygen cylinder holder;2 IV drip holes; 1 S.S. IV pole ; 1 mattress



YSHB-SE
Electric Patient Transfer Vehicle



- ★ The transfer device is used to transfer the patient automatically between the operation table (or medical bed) and the device without any move action by the operator. It can effectively avoid causing any further injury to the patient to protect the patient and reduce operator's physical intervene during the transfer course.
- ★ The device height is adjustable to any types of medical bed or operation table.
- ★ The foot brake is equipped to control 4 wheels which can be used to lock the device during the transfer procedure.
- ★ The device has retractable safety side rails on the both sides to protect the patient and it also has IV pole.
- ★ It is easy to use hand control panel to transfer the patient on the operation table or hospital bed.
- ★ Storage batteries make it safe and convenient to operate and use the transfer device.
- ★ Surface measures (L*W): 2040mm*740mm (±20mm)
- ★ Vertical movement: 570mm-870mm (±20mm)
- ★ Transfer distance of board: 510mm
- ★ Maximum load: 135KG
- ★ Storage battery: (12V, 20AH)*2



YSHB-SE-A
Electric Patient Transfer Vehicle



- ★ New concept of patient transfer-safe, comfortable, easy-operate
- ★ Special for overweight patient
- ★ The transfer device is used to transfer the patient automatically between the operation table (or medical bed) and the device without any move action by the operator. It can effectively avoid causing any further injury to the patient to protect the patient and reduce operator's physical intervene during the transfer course.
- ★ The device height is adjustable to any types of medical bed or operation table.
- ★ The foot brake is equipped to control 4 wheels which can be used to lock the device during the transfer procedure.
- ★ The device has retractable safety side rails on the both sides to protect the patient and it also has transfusion pole.
- ★ It is easy to use hand control panel to transfer the patient on the operation table or hospital bed.
- ★ Storage batteries make it safe and convenient to operate and use the transfer device.
- ★ Dimension: 2140mm*810mm
- ★ Height Range: 650mm~1050mm
- ★ Distance of Transfer Board: 600mm
- ★ Speed of transfer board: 300mm/s
- ★ Battery specification: (12V 20AH) *2PCS



Hospital Furnitures

Stretcher & Wheel Chair



YSHB-YHR-A2



- ★ The height of loading wheels is adjustable to fit different height of ambulance board.
- ★ The backrest can be adjusted from 0~ 65° angle by gas spring
- ★ Swing-down side rails to protect patient when transfer
- ★ The mechanical legs are controlled by the handles in both sides, there is a safety lever to control front and hind legs to fold while loading and prevent collapse while standing
- ★ Only one person can slide the stretcher on to the ambulance
- ★ The stretcher can be locked by the fixing device when it is in the ambulance
- ★ φ125mm wide rubber castors, two with brakes
- ★ Feature: high-low position, 5cm mattress, two straps
- ★ Material: durable aluminum alloy
- ★ Dimension: 195*55*(25-85)cm
- ★ Backrest: 0-65 degree
- ★ Loading Weight: max.160kg (350 LBS)



YSHB-YHR-A8



- ★ Used in hospital or emergency center to transport patient.
- ★ Foldable Emergency Ambulance Stretcher, made of high strength aluminum alloy,
- ★ Only one person can operate the stretcher to help patient up and down.
- ★ Folding structure will automatically react when being hit.
- ★ Feature: High-low position, 5cm mattress, two straps, I.V. Pole
- ★ Loading Weight: Max.160kg (350LBS)
- ★ Backrest: 0-88 degree
- ★ Dimension: 191x57x91cm
- ★ Folded: 206x57x28cm



YSHB-RC-B5 MRI Stretcher



- ★ Specially designed for MRI room
- ★ Adopting imported nonmagnetic castors
- ★ Product Size(L*W*H): 190*55*60cm
- ★ Self weight: 20kg
- ★ Loading Capability: 159kg
- ★ Package Quantity 1pcs
- ★ Package Size 192*64*50cm
- ★ G.W 26kg



**YSHB-YHR-LD01
Powered Stair Chair**



- ★ Material: Durable aluminum alloy
- ★ Seat: Aluminum board & spongy cushion
- ★ Loading: max.200kg (440 LBS)
- ★ Battery/Charger: 29.4V 13Ah/29.4V 2A
- ★ Motor: 24V 200W 120W
- ★ Charging Time: 6-7 hours
- ★ Dimension: 110*52*155cm
- ★ Folded: 110*52*28cm



**YSHB-YHR-W1
Stair Chair Stretcher**



- ★ This stair chair features two – 5" rear wheels that make for easy transport on most surfaces
- ★ It comes with a heavy duty vinyl protective cover that is durable and easy to clean
- ★ With 3 heavy duty quick release restraint straps to secure individuals during transport or while moving the patient up or down stairs.
- ★ Material: Durable aluminum alloy+vinyl seat
- ★ Loading weight: max.160kg (350LBS)
- ★ Dimension: 70*54*92cm
- ★ Folded: 70*54*18cm
- ★ Packing: 72*56*20cm 1 set/CTN
- ★ N.W: 9kg / G.W.: 11kg



YSHB-DJ101/102/103



- ★ YSHB-DJ101 Scoop stretcher (210*44*6cm)
- ★ YSHB-DJ102 Basket stretcher (216*61*18cm)
- ★ YSHB-DJ103 Spine board stretcher (183*45*6cm)



YSHB-WC101



- ★ Surface treatment (optional chrome-plated, spray plastic)
- ★ Seat (optional soft seat, hard seat, toilet), seat width 46cm
- ★ Small wheel 7-inch PU solid wheel, 24-inch large wheel (optional plastic body wheel, banner wheel)
- ★ Tire (optional inflatable, non-inflatable)
- ★ Load-bearing 100kg



YSHB-WC102



- ★ Adjustable leg support, surface treatment (optional chrome-plated, plastic sprayed)
- ★ Seat (optional soft seat, hard seat, toilet), seat width 46cm, small wheel 7-inch PU
- ★ Solid wheel, 24-inch large wheel (optional plastic all-in-one wheel, banner wheel), tire optional inflatable, non-inflatable)
- ★ Movable leg support
- ★ Load-bearing 100kg



**YSHB-WC46A11
Electric Wheel Chair**



- ★ Powder coating steel frame; Tube:33.5mm*15.5mm*1.5mm;
- ★ Double layers of seat cushions and backrest cushions.Upper breathable mesh fabric seat cushion with 5cm high resilient decompression sponge is detachable and washable;
- ★ 250W*2 two stage transmission variable speed helical gear motor;
- ★ Flip-up armrest with PU pad;Foldable backrest ;
- ★ Detachable legrest and foldable footplate and calf strap;
- ★ With seat belt, hand brake and anti-tippers;
- ★ Two staged shock absorbers;
- ★ Max speed: 6km/h
- ★ Range: 12±10%km
- ★ Max load: 100kg
- ★ Castor: 8"PU
- ★ Rear wheel: 12"pneumatic rear wheel
- ★ Battery: lead acid battery(24V/12AH)



Hospital Furnitures

Other Furnitures



YSHB-ZCC09/09A



- ★ Size: 1900*550*600-800mm
- ★ Backrest section can be lifted up from: 0°to 75°easily by gas spring
- ★ Foot section Adjustable by gas spring: 0-45°
- ★ Load capacity: 200KGS
- ★ Net weight: 60KGS
- ★ Gross weight: 70KGS



YSGH1101



- ★ Size (mm): 1900*600*650
- ★ Materials: #201 Stainless steel
- ★ Colors: Black/Blue
- ★ Standards: Environmentally-friendly leather/High elastic sponge
- ★ Optionals Pillow



YSGH1102-A



- ★ Size (mm): 1900*600*650
- ★ Materials: Coated steel
- ★ Colors: Black/Blue
- ★ Standards: Environmentally-friendly leather/High elastic sponge
- ★ Optionals Pillow



YSHB-KX31



- ★ Flat bedside cabinet
- ★ Size: 46.3*47.5*75cm
- ★ Material: ABS engineering plastic
- ★ Color: blue



YSHB-KX852/890



- ★ Luxury lifting dining table with pp top board (YSHB-KX852)
- ★ Size: 900x450 (height adjustment: 650x900mm)
- ★ Luxury lifting dining table with wood top board (YSHB-KX890)
- ★ Size: 800x400 (height adjustment: 650x900mm)



YSHB-JTD40



- ★ Double steps
- ★ Material: stainless steel



YSHB-PHY02



- ★ Attendant chair
- ★ Material: Coated steel with sponge
- ★ Size: 190*60*40cm
- ★ It can be extended to be a bed



YSHB-PF02/03/04



- ★ Foldable wards screen
- ★ Material: stainless steel
- ★ Size: 2-fold/3-fold/4-fold



YSHB-SYG40B



- ★ Mobile drip stand
- ★ Material: stainless steel
- ★ Size: 45*21cm



Hospital Furnitures

Other Furnitures



YSHB-HZY01



- ★ Patient waiting chair
- ★ Material: coated steel
- ★ Size: 188*68*45/80cm



YSHB-SSD02



- ★ Surgeon stool
- ★ Whole stool is made of stainless steel
- ★ With leather pad
- ★ With 5 Casters
- ★ With back rest
- ★ Gas spring control high low



YSGH1241-A



- ★ Infusion chair
- ★ Size (mm): 1500*620*400
- ★ Material: Coated steel
- ★ Color: blue
- ★ Standard: IV pole/Hand rails



YSHB-STC05



- ★ Stainless cart
- ★ Material: stainless steel
- ★ Size: 90*50*90cm



YSHB-TP10



- ★ Mobile instrument tray
- ★ Material: stainless steel
- ★ Size: 65*40*94*140cm



YSHB-ZJT01



- ★ Stainless surgical bowl
- ★ Material: stainless steel



YSHB-FP05



- ★ Stainless instrument tray
- ★ Material: stainless steel
- ★ Size: 21.5*10.5*2cm / 23*17*2.85cm / 26*20*3cm
- ★ Other sizes can be customized upon request.



YSHB-YP05



- ★ Stainless kidney tray
- ★ Material: stainless steel
- ★ Size: 20.4*13*2.5cm / 19*11.6*2.5cm / 23*12.5*2.5cm
- ★ Other sizes can be customized upon request.



YSHB-XDC01



- ★ Stainless sterilization drum
- ★ Material: stainless steel
- ★ Diameters: 17/20/23/26/30cm



YSHB-XSC02A



- ★ Stainless sink (double positions)
- ★ Material: stainless steel
- ★ Size: 160*61*80/98cm



YSHB-XSC02B



- ★ Stainless sink with mirror
- ★ Material: stainless steel
- ★ Size: 140*585*80/190cm (Other sizes customizable)



YSHB-CCG96



- ★ Stainless instrument cabinet
- ★ Material: stainless steel
- ★ Size: 96*40*175cm



Hospital Furnitures

Medical Trolley & Waste Bin



**YSHB-ET750
Emergency Trolley**



- ★ Medical emergency trolley
- ★ Material: quality engineering plastic
- ★ Size: 750*475*930mm



**YSHB-CT626
Clinical Treatment Trolley**



- ★ Clinical treatment trolley
- ★ Material: quality engineering plastic
- ★ Size: 626*475*920mm



**YSHB-AT750
Anesthesia Trolley**



- ★ Medical anesthesia trolley
- ★ Material: quality engineering plastic
- ★ Size: 750*475*930mm



YSHB-MT750 Medicine trolley



- ★ Medical medicine trolley
- ★ Material: quality engineering plastic
- ★ Size: 750*475*930mm



YSHB-NT626 Nursing Trolley



- ★ Clinical treatment trolley
- ★ Material: quality engineering plastic
- ★ Size: 626*475*930mm



**YSHB-WT480A1
Workstation Trolley**



- ★ Tabletop is made of ABS material by injection molding techniques
- ★ Height of tabletop adjustable
- ★ With space for mouse ,keyboard
- ★ Two drawers
- ★ Fan-shaped chassis, move more steady
- ★ Size: 480*520*830-1030mm



YSHB-WWC99



- ★ Laundry collecting trolley
- ★ Material: stainless steel
- ★ Size: 95*50*85



YSHB-WWC95



- ★ Laundry collecting trolley
- ★ Material: ABS
- ★ Size: 99*54*94



YSHB-STC50



- ★ Cylinder cart
- ★ Material: stainless steel
- ★ Size: 33*33*120cm



YSHB-LJT01



- ★ Medical waste bin
- ★ Material: plastic
- ★ Capacity: 15L/20L/30L/40L/50L/60L/80L/100L/120L/240L



YSHB-LQH01



- ★ Medical sharp container
- ★ Material: plastic
- ★ Capacity: 0.5L/1L/2L/3L/4L/5L/6.5L/7L/8L/10L/12L/15L/18L



Bed Pan & Urinal



- ★ Bed Pan & Urinal