

WELLMEDPRO

TREATMENT AND CARE MEDICAL PRODUCTS



Nebulizer | Nebulizer Cup | Nasal Aspirator | Nasal Irrigator

Breast Pump | Insect Bite Healer | Otoscope | TENS&EMS | Ultrasonic Cleaner



About Wellmedpro

Your one-stop equipment provider

Wellmedpro provides medical and veterinary equipment, with headquarters in China and a Latin America sales office in Panama.

Sourcing and commercial support

Contact our team for product selection, quotations, availability and coordination of orders for your market.



Product documentation

Documentation for your selected products

Ask Wellmedpro for technical specifications and available regulatory documentation for the exact model and destination market.

Project requirements

Share your required quantities, intended use and delivery destination so our team can prepare a tailored proposal.



From selection to delivery

1. Define your requirements

Choose products and share quantities, technical requirements and delivery destination.

2. Request a proposal

Our team coordinates pricing, availability and product documentation.

3. Coordinate your order

Contact Wellmedpro to review order details and delivery arrangements.

Contents

Focusing On Home Use Medical Products

.....

Nebulizer	01-09
-----------	-------

.....

Nebulizer Cup	10-11
---------------	-------

.....

Breast Pump	12-21
-------------	-------

.....

Insect Bite Healer	22-28
--------------------	-------

.....

Nasal Aspirator	29-32
-----------------	-------

.....

Nasal Irrigator	33-36
-----------------	-------

.....

Otoscope	37-39
----------	-------

.....

TENS&EMS	40-41
----------	-------

.....

Ultrasonic Cleaner	42
--------------------	----

NEBULIZER



Compressor Nebulizer

Efficient Delivery to Target Site

High-Performance Silent Compressor



NEB-CP01

Compressor Nebulizer

- Set the time 0-30 minutes / Always on
- Powerful compressor, thorough nebulization
- Adjustable nebulization rate

Particle size(MMAD): $2.2 \mu\text{m} \pm 10\%$

Nebulization Rate: $\geq 0.3 \text{ mL/min}$

The proportion of aerosol particles smaller than $\Phi 5\mu\text{m}$ is $\geq 81.24\%$

Size: 194 x 179 x 108 mm

Power source: 230V~50 Hz

Noise level: $\leq 53 \text{ dB}$ (No liquid)

Weight: approx. 1.5kg



Anti-slip
foot mats



Adjustable particles



For children
& adults



Heat exhaust port



Replacement
Air Filters



Timing function &
Always on



Quiet motor

Two Modes of flow settings



NEB-CP02

Compressor Nebulizer

- Two Modes: Tender for sensitivity, Power for depth
- Higher efficiency for superior atomization
- Internal dual shock absorption for effective noise reduction
- For a variety of professional and home settings

Particle size(MMAD): $2.8 \mu\text{m} \pm 10\%$ (Grade 1) , $2.2 \mu\text{m} \pm 10\%$ (Grade 2)

Nebulization Rate: $\geq 0.34 \text{ mL/min}$

Grade 1: Proportion of aerosol particles smaller than $\Phi 5\mu\text{m} \geq 75.50\%$

Grade 2: Proportion of aerosol particles smaller than $\Phi 5\mu\text{m} \geq 80.90\%$

Size: 194 x 179 x 108 mm

Power source: 230V~50 Hz

Noise level: $\leq 53 \text{ dB}$ (No liquid)

Weight: approx.1.5kg



Adjustable Rate



Adjustable particles



For children
& adults



Anti-slip
foot mats



Replacement
Air Filters



Timing function &
Always on



Heat exhaust port



Quiet motor

Customized Timing for Targeted Therapy



NEB-CP03

Compressor Nebulizer

- Set the time 0-30 minutes / Always on
- Micron atomization for deep lung penetration
- Scientific grille enhances thermal stability

Particle size(MMAD): $2.2 \mu\text{m} \pm 10\%$

Nebulization Rate: $\geq 0.3 \text{ mL/min}$

The proportion of aerosol particles smaller than $\phi 5\mu\text{m}$ is $\geq 80.28\%$

Size: 176 x 159 x 98 mm

Power source: 230V~50 Hz

Noise level: $\leq 53 \text{ dB}$ (No liquid)

Weight: approx. 1.3kg



Anti-slip
foot mats



Adjustable particles



For children
& adults



Heat exhaust port



Replacement
Air Filters



Timing function &
Always on



Quiet motor

Two Modes of flow settings



NEB-CP04

Compressor Nebulizer

- Adjustable droplet size, precise treatment
- Powerful compressor, thorough nebulization
- Adjustable nebulization rate
- Silencer chamber designed with fluid dynamics

Particle size(MMAD): $2.8 \mu\text{m} \pm 10\%$ (Grade 1), $2.2 \mu\text{m} \pm 10\%$ (Grade 2)

Nebulization Rate: $\geq 0.34 \text{ mL/min}$

Grade 1: Proportion of aerosol particles smaller than $\Phi 5\mu\text{m} \geq 75.50\%$

Grade 2: Proportion of aerosol particles smaller than $\Phi 5\mu\text{m} \geq 80.90\%$

Size: 178 x 159 x 97 mm

Power source: 230V~50 Hz

Noise level: $\leq 53 \text{ dB}$ (No liquid)

Weight: approx.1.3kg



Adjustable Rate



Adjustable particles



For children
& adults



Anti-slip
foot mats



Replacement
Air Filters



Timing function&
Always on



Heat exhaust port



Quiet motor

NEB-CP05

Compressor Nebulizer

- One-Touch Continuous Mode: For sustained atomization
- Higher efficiency for superior atomization
- Internal dual shock absorption for effective noise reduction
- For a variety of professional and home settings

Particle size(MMAD): $2.2 \mu\text{m} \pm 10\%$

Nebulization Rate: $\geq 0.34 \text{ mL/min}$

The proportion of aerosol particles smaller than $\phi 5\mu\text{m}$ is $\geq 80.28\%$

Size: 178 x 159 x 97 mm

Power source: 230V~50 Hz

Noise level: $\leq 53 \text{ dB}$ (No liquid)

Weight: approx. 1.4kg



Adjustable Rate



Adjustable particles



For children
& adults



Anti-slip
foot mats



Replacement
Air Filters



Timing function &
Always on



Heat exhaust port



Quiet motor



NEB-MP01

Compressor Nebulizer

- Micron atomization, high respirable particle ratio
- Reliable, portable, and stylish cartoon design
- Comfortable aerosol, no coughing sensation

Particle size(MMAD): $3.3 \mu\text{m} \pm 10\%$

The proportion of aerosol particles smaller than $\phi 5\mu\text{m}$ is $\geq 40\%$

Size: 136 x 103.6 x 68.3 mm

Power: 10 W

Noise level: $\leq 58 \text{ dB}$ (at 1m distance)

Weight: approx. 255g



Adjustable Rate



Heat exhaust port



For children
& adults



Anti-slip
foot mats



Replacement
Air Filters

Fine Particles & Efficient Treatments



NEB-MP02

Compressor Nebulizer

- Micron atomization, high respirable particle ratio
- Robust, low-maintenance, long working life
- Comfortable aerosol, no coughing sensation

Particle size(MMAD): $3.3 \mu\text{m} \pm 10\%$

The proportion of aerosol particles smaller than $\phi 5\mu\text{m}$ is $\geq 40\%$

Size: 130 x 132.7 x 75.5 mm

Power: 10 W

Noise level: $\leq 58 \text{ dB}$ (at 1m distance)

Weight: approx. 310g



Adjustable Rate



Heat exhaust port



For children
& adults



Anti-slip
foot mats



Replacement
Air Filters

Compact Design, Convenient to carry



NEB-MP03

Compressor Nebulizer

- Micron atomization, high respirable particle ratio
- Reliable, portable, and stylish cartoon design
- Comfortable aerosol, no coughing sensation

Particle size(MMAD): $3.3 \mu\text{m} \pm 10\%$

Nebulization Rate: $\geq 0.25 \text{ mL/min}$

The proportion of aerosol particles smaller than $\phi 5\mu\text{m}$ is $\geq 40\%$

Size: 124 x 70 x 53 mm

Power: 10 W

Noise level: $\leq 58 \text{ dB}$ (at 1m distance)

Weight: approx. 200g



Adjustable Rate



Heat exhaust port



For children
& adults



Anti-slip
foot mats



Replacement
Air Filters



NEBULIZER CUP

- 3-in-1 kit: Child mask, adult mask & mouthpiece for all ages.
- One-handed 3-stage mist control for adjustable output.
- Particle size(MMAD) 2.5 μ m mist ($\pm$ 10%) for gentle inhalation.
- $\leq$ 0.6mL residual liquid for efficient medication use.
- >80% of particles under 5 μ m for effective delivery.
- Fully detachable for easy cleaning and home reuse.
- Medical-grade & phthalate-free; passed safety tests.
- Seamless core with no adhesives for reliable performance.
- Size: 94 × 39 × 95 mm

NEK-01

- High Pressure, Extended Lifespan
- Service Life: 1 year or 300 cycles (whichever comes first).
For use with AC atomizers.



Compatible Accessories

AeroJet-01

- Low Pressure, Limited Lifespan
- Service life: 1 year or 300 cycles. Compatible with DC atomizers



NEK-Mini-01

- High Pressure, Limited Lifespan
- Service Life: 15 cycles. For use with AC atomizers.




AeroJet-Mini-01

- Low Pressure, Extended Lifespan
- Service life: 15 cycles. For use with DC atomizers



BREAST PUMP



Wearable Breast Pump



BRP-301

Electric Breast Pump

- Proprietary pump tech with multiple patents
- Max setting, lasts 1.8 hours
- 2 modes, 9 adjustable vacuum levels

LED screen | 6oz/180mL capacity

Power Source: D.C. 3.7V/Li-battery 1000mAh

Size: 110 x 73.8 x 130.3 mm (4.3 x 2.9 x 5.1 inches)

Screen Size: 35 x 18.6 mm (1.56 inches)

Weight: About 225g

Battery Life: 1.8h (Max setting)



Designed for
Multitasking



Super large
LED screen



Massage &
Expression



Time display



2 Modes
& 9 Levels



Food-grade
silicone

Leak-proof Collection Cup Design



BRP-302

Electric Breast Pump

- Upgraded diaphragm for strong, lasting suction
- Soft touch, less friction on the skin
- Two working modes, 9 adjustable levels

LED screen | 6oz/180mL capacity

Power Source: D.C.3.7V/Li-battery 1000mAh

Size: 65.7 x 108.6 x 125.6 mm (2.5 x 4.3 x 4.9 inches)

Screen Size: 27 x 18 mm (1.28 inches)

Weight: About 215g

Battery Life: 2h (Max setting)



Anatomically
design



2 Modes
& 9 Levels



Food-grade
silicone



Night light



Anti-backflow design



Ultra-light
No Burden

Wearable Breast Pump



BRP-306

Electric Breast Pump

- Proprietary pump tech with multiple patents
- High-frequency lactation simulation
- Vertical placement to avoid tipping

LED screen | 6oz/180mL capacity

Power Source: D.C.3.7V/Li-battery 1000mAh

Size: 74 x 111 x 141mm (2.9 x 4.3 x 5.5 inches)

Screen Size: 27 x 18mm (1.28 inches)

Weight: About 195g

Battery Life: 1.8h(Max setting)



Designed for
Multitasking



Super large
LED screen



Massage &
Expression



Time display



2 Modes
& 9 Levels



Food-grade
silicone

BRP-601

Dual-side Wearable Breast Pump

- The breast cups can be controlled independently
- The design of a wearable hands-free electric breast pump
- Easy to clean and fully detachable
- The run time can be set



LED screen | 6oz/180mL capacity

Power Source: Li-battery 2000mAh

Size: 160 x 90 x 42 mm (6.3 x 3.54 x 1.65 inches)

Screen Size: 43.5 x 77 mm (3.48 inches)

Weight: About 500g

Battery Life: 2h (Max setting)



Ultra-light
No Burden



Intuitive control panel



High capacity
Li-Battery



Setting the
runtime



Designed for
flexibility



Double
pumping

BREAST PUMP SET

SET 01

3 in 1 set



BRP-101

BRP-304

BRP-301

SET 02

3 in 1 set



BRP-101

BRP-305

BRP-306

SET 03

2 in 1 set



BRP-304

BRP-301

SET 04

2 in 1 set



BRP-305

BRP-306

IF Design Award



BRP-201

Dual-side Wearable Breast Pump

- Mimics baby's natural sucking rhythm
- Efficient silent pump for quick response
- More milk, lightweight and comfortable
- Smart dual-side remote control

LED screen | 6oz/180mL capacity

Power Source: Li-battery 2000mAh

Size: 170 x 80 x 41.7 mm (3.6 x 2.6 x 1.2 inches)

Screen Size: 43.5 x 77 mm (3.48 inches)

Weight: About 350g

Battery Life: 2.5h(Max setting)



Ultra-light
No Burden



Intuitive control panel



High capacity
Li-Battery



Setting the
runtime



Designed for
flexibility



Double
pumping

BRP-202

Dual-side Wearable Breast Pump

- The flexibility to pump any time, anywhere
 - Designed for multi-tasking
 - Fits most nursing bras
- Vertical placement to avoid tipping



BRP-203

Dual-side Breast Pump

- 4 parts, easy to clean and maintain
- Traditional bottle collection kit
- Customized waveform, e.g. automatic, mixed mode, etc



Connecting Disposable Milk Bags



Storage in one step



Clear dual-color alert

BRP-303

Electric Breast Pump

- Advanced pump tech with multiple patents
- Instant storage with disposable milk bags
- Precise negative pressure control

LED screen | 5oz/150mL capacity

Power Source: Li-battery 1100mAh

Size: 126 x 106 x 67 mm (4.96 x 4.17 x 2.64 inches)

Battery Life: 1.8h(Max setting)



Ultra-light
No Burden



Micro-vibration tech



Food-grade
silicone




Dual-color LED
mode indicator



3 Modes
& 9 Levels



Anti-backflow Design



Quick Treatment for Insect Bite

IBR-10 / IBR-20

Insect Bite Heater

- High-precision temperature sensor
- Professional algorithm controls
- Two settings for child & adult

Measurement Range: 47°C/116.6°F, 50°C/122°F

Measurement Accuracy: $\pm 2^{\circ}\text{C}$

Power Source: 2 x batteries, type AAA

Size: 109 x 34 x 25 mm (4.29 x 1.34 x 0.98 inches)

Battery Life: 340 times (First gear)

Weight: About 25g



Cater to
different skin



MD certification
cleared



Auto power off



Fast-acting
relief



Locked mode
(optional)



Compact & Portable

IBR-30 / IBR-40

Insect Bite Healer

- Two settings for child & adult
- Professional algorithm controls
- Safe and non-invasive relief

Measurement Range: 47°C/116.6°F, 50°C/122°F

Measurement Accuracy: $\pm 2^{\circ}\text{C}$

Power Source: 2 x 1.5V AAA batteries

Size: 109 x 34 x 25 mm (4.29 x 1.34 x 0.98 inches)



Compact & Portable



MD certification
cleared



Auto power off



Fast-acting relief



Locked mode
(optional)



Cater to different skin



One-Click Operation



Compact & Portable



Highly waterproof

IBR-50 / IBR-60

Insect Bite Healer

- Highly Waterproof, outdoor-friendly
- Ergonomic design for optimal contact
- Professional algorithm controls

Measurement Range: 47°C/116.6°F, 50°C/122°F

Measurement Accuracy: $\pm 2^{\circ}\text{C}$

Power Source: 2 x 1.5V AAA batteries

Size: 135 x 23 x 23 mm (5.31 x 0.91 x 0.91 inches)



Pocket-sized



Highly waterproof



Auto power off



Fast-acting relief



Locked mode
(optional)



Cater to different skin



Patent No.
ZL 2023 1 1636339.6
ZL 2023 2 1799078.5



IBR-70 / IBR-75 / IBR-80

Insect Bite Heater

- **Ultra-Portable:** Fits anywhere, weighs nothing
- **Universal Use:** Plugs into phones, power banks, laptops
- **Dual Heat Modes:** Gentle 3S or Powerful 6S
- **Smart & Safe:** Precise temperature control with overheating protection

Display Type: Tri-color Indicator Light

Auto Standby Time: 1 minute

Power Source: Phone's Type-C Port

Size: 15.5 x 30 x 7.5 mm (0.61 x 1.18 x 0.30 inches)



Pocket-sized



Auto power off



Fast-acting relief



Auto Standby Time
1 minute



Cater to different skin



IBR-600

Rechargeable battery Insect Bite Heater

- **Compact & Portable:** Pocket-sized for on-the-go ease
- **Dual Heating Modes:** Gentle 3S for sensitivity, Powerful 6S for deep relief
- **Precise Temp Control:** Smart sensor for safety and consistent warmth

Power Source: DC 3.7V Li-battery 200mAh

Size: 18.5 x 77 mm (0.73 x 3.03 inches)

Weight: About 16.8g



Rechargeable Battery



IP65 Water



Pocket-sized



Auto power off



Fast-acting relief



Cater to different skin



Sleep mode

IBR-800

Rechargeable battery Insect Bite Healer



- **Compact & Portable:** Pocket-sized for on-the-go ease
- **Dual Heating Modes:** Gentle 3S for sensitivity, Powerful 6S for deep relief
- **Precise Temp Control:** Smart sensor for safety and consistent warmth

Power Source: DC 3.7V Li-battery 300mAh

Size: 40.3 x 30.5 x 52 mm (1.59 x 1.20 x 2.05 inches)

Weight: About 23.5g



Rechargeable Battery



Pocket-sized



Sleep mode



Auto power off




Fast-acting relief



Cater to different skin

NASAL ASPIRATOR



Reliable and Adjustable



Four nose tip sizes



ASP-401

Nasal Aspirator

- Precise negative pressure control
- One-touch operation
- Ergonomic design, Any-angle use

Power Source: 2 x1.5V AA batteries

Suction power : 26kPa & 46kPa

Size: 155 x 56 x 44 mm (6.1 x 2.2 x 1.73 inches)



Anti-reflux, no
reverse suction



Easy to disassemble
& clean



Food-grade
Silicone



Dual suction
levels



Self-standing base



One-touch operation

Controllable Suction Levels



Four nose tip sizes



Color lights

ASP-402

Nasal Aspirator

- Five-Level suction adjustment
- Three soothing music options
- Attention-grabbing color lights

Power Source: 3.7V Li-battery 450mAh

Suction power: -55~-75kPa

Size: 168 x 54 x 47 mm (6.61 x 2.13 x 1.85 inches)



Rechargeable



Easy to
disassemble & clean



Food-grade
silicone



Large LED
display



Anti-reflux, no
reverse suction



One-touch operation



ASP-101

Nasal Aspirator Accessory

- Separated from the main unit with low noise.
- Fully washable for easy cleaning.
- Flexible combination: one main unit with two accessories for dual functions (nasal cleansing and suction).
- Compact and user-friendly, suitable for use while standing, lying down, or sitting in various scenarios.



low noise



Easy-to-clean



Dual Function



User-friendly Design

NASAL IRRIGATOR





ANI-801 / ANI-802

Nasal Irrigator

- Gentle spray, easy on the nose
- Simple one-button operation for all ages
- Choose between pulse or continuous flow modes
- Separate tanks for clean solution and waste

Air Pump (Big): DC 6V, Maximum inflation pressure $\geq 93\text{kPa}$

Air Pump (Small) : DC 9V, Maximum inflation pressure $\geq 55\text{kPa}$

Power Source: 3.7V Li-battery 800mAh

Spray Rate (Big): Continuous Mode $\geq 15.8\text{ ml/min}$, Pulse Mode $\geq 6.3\text{ ml/min}$

Spray Rate (Small): Continuous Mode $\geq 14.4\text{ ml/min}$, Pulse Mode $\geq 5.3\text{ ml/min}$

Size: 64 x 52 x 167 mm (2.52 x 2.05 x 6.58 inches)



Anti-reflux, no
reverse suction



Easy to disassemble
& clean



Food-grade
Silicone



Self-standing base



One-touch operation



Type-C charging

ANI-803 / ANI-804

Nasal Irrigator

- Gentle spray, easy on the nose
- Simple one-button operation for all ages
- Choose between pulse or continuous flow modes
- Separate tanks for clean solution and waste

Air Pump (Big): DC 6V, Maximum inflation pressure $\geq 93\text{kPa}$

Air Pump (Small) : DC 9V, Maximum inflation pressure $\geq 55\text{kPa}$

Power Source: Powered by DC 5V via plug

Spray Rate (Big): Continuous Mode $\geq 15.8\text{ ml/min}$, Pulse Mode $\geq 6.3\text{ ml/min}$

Spray Rate (Small): Continuous Mode $\geq 14.4\text{ ml/min}$, Pulse Mode $\geq 5.3\text{ ml/min}$

Size: 64 x 52 x 167 mm (2.52 x 2.05 x 6.58 inches)



Anti-reflux, no reverse suction



Easy to disassemble & clean



Food-grade Silicone



Self-standing base



One-touch operation



ANI-101

Nasal Irrigator Accessory

- Separated from the main unit with low noise.
- Fully washable for easy cleaning.
- Flexible combination: one main unit with two accessories for dual functions (nasal cleansing and suction).
- Compact and user-friendly, suitable for use while standing, lying down, or sitting in various scenarios.



low noise




Easy-to-clean



Dual Function



User-friendly Design



4-in-1 Multifunctional Otoscope



Disposable Probe Cover



Inspection for Ear,
Nose, Throat & Mouth

OT-501

Otoscope

- Ultra-clear real-time video
- Waterproof lens
- Intelligent visible APP

Battery Capacity: 800mAh Li-battery

Input Current: USB 5V/1A

Size: 23 x 55 x 155 mm (0.9 x 2.2 x 6 inches)



3MP HD
camera



3-Axis Gyroscope



High-speed WIFI chip



4mm
Ultra-slim lens



Cold 6 LED Lights



360° Panorama

Lightweight and Portable



EWR-10 / EWR-50

Intelligent Visual Ear Wax Removal

- Ultra-clear real-time video
- Food grade soft silicone tip
- Intelligent visible APP

Battery Capacity: 3.7V 200mAh Li-battery

Input Current: DC 5V/1A

Size: 14 x 14 x 139 mm (0.55 x 0.55 x 5.47 inches)



3MP HD camera



3-Axis Gyroscope



High-speed WIFI chip




3.5mm Ultra-slim lens



Cold 6 LED Lights



360° Panorama





EPT-601 / EPT-602

Low frequency electronic pulse Therapeutic Device

- Rapid Pain Relief: Targets pain at its source.
- Improved Circulation: Boosts blood flow and nutrient delivery
- Restored Function: Eases stiffness and muscle atrophy
- Reduced Inflammation: Effectively diminishes swelling and pain
- Accelerated Recovery: Shortens healing time

Battery Capacity: 3.7V 200mAh Li-battery

Charging Voltage: DC5V

Rated Power: 3.7W

Net Weight: 120g

Size: 132 x 64 x 14mm (5.20 x 2.52 x 0.55 inches)

Patch Size: 40 x 40mm

Output Cable Length: 1.5m

Operating Temperature: 5°C to 40°C



3 Massage Modes



For 6 Body Areas



Rechargeable & Long-lasting



Washable
Electrode Pads



18 Adjustment Options

USC-10

Ultrasonic Cleaner

- Large-Capacity SUS304 Cleaning Tank
- Viewing window for monitoring the ultrasonic cleaning process
- Automatic shut-off for eye safety, with intuitive touch-button status indicators



Size: 101 × 116 × 86 mm

Housing Material: ABS + PC

Tank Material: SUS304 Stainless Steel

Net Weight: 310 g

Tank Capacity: 200 ml

Ultrasonic Transducer Size: Ø 50 × 2.6 mm

Ultrasonic Transducer Power: 20W ± 3W

Ultrasonic Frequency: 45 kHz ± 3 kHz

UVC Lamp Power: 0.5 W



Ultrasonic Frequency
45 kHz



Tank Material
SUS304 Stainless Steel



Tank Capacity 200 ml



Ultrasonic Transducer
Power 20W




Contact Wellmedpro

Website: www.wellmedpro.com

Panama / WhatsApp: +507 6351 5786

China headquarters: +86 136 3275 4074

LinkedIn: [linkedin.com/company/wellmed-pro/](https://www.linkedin.com/company/wellmed-pro/)



WhatsApp