




Wellmedpro

Room 1617-18 Star House 3 Salisbury Road Tsim Sha Tsui Kowloon Hong
KongTel.:+86-755-82830455 / +86-755-82041897
www.wellmedpro.com

WELLMEDPRO



Wellview Universal Ltd

Contents

Structure	1	▶
Electric Wheelchair	3	▶
Cerebral Palsy Wheelchair	5	▶
Leisure Type Sports Wheelchair	6	▶
Aluminum transpot wheelchair	7	▶
Aluminum transport lightweight	8	▶
Functional Child Wheelchair	11	▶
Care aluminum wheelchair	13	▶
Aluminum Manual Wheelchair	15	▶
Functional Aluminum Manual Wheelchair	17	▶
Steel Manual Wheelchair	21	▶
Functional Steel Manual Wheelchair	26	▶
Functional Commode Wheelchair	31	▶
Commode Wheelchair	36	▶
Shower Chair	41	▶
Walker	42	▶
Walking stick	45	▶
Accessories	47	▶
Fabric	50	▶

Wheelchair



WHEEL CHAIR MODEL FUNCTION CODE

	Remark
A	DROP BACK HANDLE
B	Mag WHEEL
C	ELEVATING LEGREST
G	RECLINING HIGH BACK
J	UNITED BRAKE
L	ALUMINUM
S	STAINLESS STEEL
Y	HARD SEAT CUSHION
Q	QUICK RELEASE
P	PU TYRE

For: ALK 908 L A B J-46

SEAT WIDTH 46cm
UNITED BRAKE
Mag WHEEL
DROP BACK HANDLE
ALUMINUM

(-35) : SEAT WIDTH 35cm (14")
(-41) : SEAT WIDTH 35cm (16")
(-43) : SEAT WIDTH 35cm (17")
(-46) : SEAT WIDTH 35cm (18")
(-51) : SEAT WIDTH 35cm (20")
(-56) : SEAT WIDTH 35cm (22")



Model / **ALK158**
Name / Electric Wheelchair

- 1) Spray power frame with breathable mesh cushion
- 2) 8" big castor, 12" PU rear wheel
- 3) Detachable plate, the armrest can be turned backwards
- 4) Can go up and down and walk upright by electric motor
- 5) With safe-light

Model	Size (cm)	→←	↔	↕	↕	↕	↕	↕	↕	↕	↕	↕	↕	↕
ALK158	61	-	44	32	20	57	105	115	43	48	110	105x61.5x65	80	98
ALK140LA	58	43	46	50	20	52	100	103	40	46	110	76x44x79	52	54
ALK110A	57/62	35	41/46	32	20	49	90	100	40	40	110	75x38x72	57	59
ALK111A	62	35	46	32	20	53	99	98	40	42	110	82x39x89	53	55
ALK122A	66	28	46	32	20	54	100	99	40	45	110	78x28x92	44	46

Model / **ALK140LA**
Name / Electric Wheelchair

- 1) Aluminum lightweight frame with comfortable cushion
- 2) 8" big castor, 20" PU rear wheel
- 3) Detachable footrest, folding backrest, the armrest can be turned backwards
- 4) Can walk by electric motor



Model / **ALK111A**
Name / Electric Wheelchair

- 1) Spray power frame with comfortable cushion
- 2) 8" big castor, 12" PU rear wheel
- 3) Detachable footrest, folding backrest, the armrest can be turned backwards
- 4) Can walk by electric motor



Model / **ALK110A**
Name / Electric Wheelchair

- 1) Spray power frame with comfortable cushion
- 2) 8" big castor, 12" PU rear wheel
- 3) Detachable footrest, folding backrest, the armrest can be turned backwards
- 4) Can walk by electric motor



Model / **ALK122A**
Name / Electric Wheelchair

- 1) Powder coated steel folding frame
- 2) Flip up armrest, detachable footrest
- 3) 200*50PU front castor, 12" PU rear wheel





Model / **ALK958LC**

Name / Cerebral Palsy Wheelchair

- 1) Aluminum lightweight frame with comfortable seat, reclining high back, adjustable headrest
- 2) 6" PVC castor, 16" wheel
- 3) Adjustable and detachable armrest, elevating footrest, seat and backrest adjustable
- 4) With safety belt



Model / **ALK959L**

Name / Cerebral Palsy Wheelchair

- 1) Powder coated aluminum lightweight frame
- 2) Adjustable and detachable armrest
- 3) Fixed footrest and adjustable footplate
- 4) 6"PVC front castor,16" PU rear wheel



Model / **ALK278L**

Name / Leisure Type Sports Wheelchair

- 1) Powder coated aluminum lightweight frame
- 2) Fixed footrest
- 3) 6"PVC front wheel, 24"pneumatic rear wheel



Model / **ALK221LP**

Name / Leisure Type Sports Wheelchair

- 1) Folding aluminum frame with fixed armrest
- 2) 3" solid front castor and 24" PU rear wheel

Model	Size (cm)	Seat Width	Seat Depth	Seat Height	Wheel Width	Wheel Diameter	Wheelbase	Wheel Offset	Armrest Height	Backrest Height	Weight (kg)	Net Weight (kg)	Gross Weight (kg)
ALK958LC	56/51/48	26	46/41/36	38	15	51	93	103	36	47	75	77x28x92	23 25
ALK959L	56/51/48	26	46/41/36	38	15	51	93	103	36	47	75	85x28x92	15 17

Model \ Size (cm)													
ALK278L	63		36	59	12.5	51	75	77	36	30	100	69x65x49	14 16
ALK221LP	59	28	36	60	7.5	45	73	90	38	31	100	82x31x74	10 12



Model / **ALK920L**

Name / Transport Wheelchair

- 1) Powder coated aluminum lightweight frame
- 2) Flip up armrest, folding footrest
- 3) 5"PVC front castor, 8"PVC rear wheel



Model / **ALK901LAJ**

Name / Transport Wheelchair

- 1) Nano Dryplating aluminum frame with brake
- 2) Drop back handle, tip-back armrest and footrest
- 3) 5" solid castor and 8" rear wheel



Model / **ALK900L**

Name / Transport Wheelchair

- 1) Collapsible portable aluminum frame,with brake
- 2) Tip-back armrest and footrest
- 3) 5" solid castor and 8" rear wheel,safety belt



Model / **ALK863LAB**

Name / Transport Wheelchair

- 1) Aluminum lightweight frame
- 2) Drop back handle, fixed armrest and tip-back footrest
- 3) 6" solid castor and 12" PU rear wheel



Model / **ALK863LABJ**

Name / Transport Wheelchair

- 1) Collapsible portable aluminum frame,with united brake
- 2) Drop back handle, fixed armrest and tip-back footrest
- 3) 6" solid castor and 12" PU rear wheel



Model	Size (cm)	Seat Width	Seat Depth	Backrest Width	Backrest Height	Wheel Width	Wheel Depth	Wheel Diameter	Wheel Weight	Weight	Capacity	Dimensions (cm)	NW	GW
ALK920L	55	33	35	20	12.5	49	91	82	32	38	75	55x33x75	6	8
ALK900L	46	33	30	20	12.5	44	91	85	33	38	75	46x33x71	7	8

Model	Size (cm)	Seat Width	Seat Depth	Backrest Width	Backrest Height	Wheel Width	Wheel Depth	Wheel Diameter	Wheel Weight	Weight	Capacity	Dimensions (cm)	NW	GW
ALK901LAJ	43/48	21	36/41	20	12.5	45	86	95	36	38	75	77x21x65	6	8
ALK863LAB	56/61	27	41/46	30	15	46	84	99	36	42	75	73x26x71	9	11
ALK863LABJ	56/61	27	41/46	30	15	46	84	99	36	42	75	73x26x71	9	11



Model / **ALK863LAJP-20"**
Name / Transport Wheelchair

- 1) Nano Dryplating aluminum frame with united brake
- 2) Drop back handle, fixed armrest and tip-back footrest
- 3) 6" solid castor and 20"PU rear wheel with aluminum frame

Model	Size (cm)														
ALK863LAJP-20"	57/62	27	41/46	50	15	46	88	100	36	42	75	81x28x71	10	13	



Model / **ALK902LABJ**
Name / Transport Wheelchair

- 1) Nano Dryplating aluminum frame with brake
- 2) Drop back handle, tip-back armrest and footrest
- 3) 5" solid castor and 12"PU rear wheel



Model / **ALK801LJP**
Name / Functional Child Wheelchair

- 1) Nano Dryplating aluminum frame with united brake
- 2) Fixed armrest and footrest
- 3) 6" solid castor and 20" PU rear wheel with aluminum frame

Model	Size (cm)														
ALK902LABJ	45/50	24	36/41	30	12.5	45	85	95	36	38	75	70x24x68	10	12	
ALK801LJP	52/57	25	36/41	50	15	49	86	86	36	39	100	76x25x88	10	12	



Model / ALK907LAP-36

Name / Functional Child Wheelchair

- 1) Nano Dryplating aluminum frame
- 2) Drop back handle, tip-back armrest and fixed footrest
- 3) 6" solid castor and 22" PU rear wheel with aluminum frame

Model / ALK802-35

Name / Functional Child Wheelchair

- 1) Nano Dryplating steel frame
- 2) Fixed armrest and footrest
- 3) 6" solid front castor, 22" solid rear wheel



Model / ALK801LAJP

Name / Functional Child Wheelchair

- 1) Aluminum lightweight frame
- 2) Fixed armrest and footrest
- 3) 6" solid castor and 20" PU rear wheel with aluminum frame
- 4) Drop back handle with united brake



Model / ALK801LAJ-12"

Name / Functional Child Wheelchair

- 1) Aluminum lightweight frame
- 2) Drop back handle with united brake, fixed armrest and footrest
- 3) 6" solid castor and 12" PU rear wheel



Model / ALK865LABJ

Name / Functional Child Wheelchair

- 1) Collapsible portable aluminum frame
- 2) Drop back handle with united brake, fixed armrest and fixed footrest
- 3) 6" solid castor and 12" PU rear wheel

Model	Size (cm)													
ALK907LAP-36	52	27	36	56	15	50	89	99	36	43	100	85x27x70	12	14
ALK907LAP-46	62	28	46	59	15	51	90	99	40	40	100	88x28x70	13	15
ALK802-35	54	23	35	56	15	47	85	97	36	42	75	87x22x85	14	16

Model	Size (cm)													
ALK801LAJP	52/57	25	36/41	50	15	49	88	86	36	39	100	76x25x71	11	13
ALK801LAJ-12"	47/52	22	36/41	30	15	49	89	85	36	39	75	68x22x71	11	13
ALK865LABJ	52/57	22	41/46	30	15	45	88	95	40	40	100	80x22x71	11	13



Model / **ALK867LB**

Name / Care Wheelchair

- 1) Nano Dryplating Aluminum lightweight frame; double seat and back (net seat + detachable soft cushion) inner net seat for summer use, outer detachable cushion for winter use
- 2) Fixed armrest and footrest
- 3) 6" PVC castor, 16"PU rear wheel



Model / **ALK972LBJ**

Name / Care Wheelchair

- 1) Nano Dryplating,Aluminum lightweight frame,with united brake
- 2) Fixed armrest and detachable footrest
- 3) Solid 8" PVC castor and 16" PU rear wheel



Model / **ALK867LABJ**

Name / Care Wheelchair

- 1) Aluminum lightweight frame;double seat and back (net seat + detachable soft cushion) inner net seat for summer use, outer detachable cushion for winter use
- 2) Drop back handle with united brake, fixed armrest and footrest
- 3) 6" PVC castor,16"PU rear wheel



Model / **ALK867LBJ**

Name / Care Wheelchair

- 1) Nano Dryplating,Aluminum lightweight frame with united brake, double seat and back(net seat + detachable soft cushion) inner net seat for summer use, outer detachable cushion for winter use
- 2) Fixed armrest and footrest
- 3) 6" PVC castor,16"PU rear wheel



Model / **ALK973LABJ**

Name / Care Wheelchair

- 1) Nano Dryplating,aluminum lightweight frame
- 2) Drop back handle with united brake
- 3) Tip-back armrest and detachable footrest
- 4) Solid 6" PVC castor and 16"PU rear wheel



Model	Size (cm)														
															
ALK867LB	51/56	23	41/46	38	15	49	85	92	40	42	100	94x23x87	10	12	
ALK867LABJ	51/56	23	41/46	38	15	49	88	94	40	42	100	88x23x73	10	12	
ALK867LBJ	51/56	23	41/46	38	15	49	85	92	40	42	100	94x23x87	10	12	

Model	Size (cm)											
-------	-----------	--	--	--	--	--	--	--	--	--	--	--



Model / ALK864LJP

Name / Aluminum Manual Wheelchair

- 1) Nano Dryplating aluminum lightweight frame with united brake
- 2) 8" solid castor and 24"PU rear wheel
- 3) Fixed armrest and footrest

Model / **ALK864LB**

Name / Aluminum Manual Wheelchair

- 1) Nano Dryplating aluminum lightweight frame
- 2) Fixed armrest and footrest
- 3) 8" solid castor and 24"PU mag rear wheel



Model \ Size (cm)														
ALK864LJP	61/66	28	41/46	59	20	46	87	103	40	43	100	92x28x90	13	15
ALK864LB	61/66	28	41/46	60	20	46	87	103	40	43	100	92x28x90	13	15
ALK864LBJ	61/66	28	41/46	60	20	46	87	103	40	43	100	92x28x90	13	15



Model / **ALK868LAJP**

Name / Aluminum Manual Wheelchair

- 1) Nano Drying, Aluminum lightweight frame; double seat and back (net seat + detachable soft cushion) inner net seat for summer use, outer detachable cushion for winter use
- 2) Drop back handle with united brake, fixed armrest and footrest
- 3) 8" PVC castor, 24" PU rear wheel with aluminum frame



Model / ALK868LP

Name / Aluminum Manual Wheelchair

- 1) Nano Dryplating, Aluminum lightweight frame
- 2) Fixed armrest and footrest
- 3) 8" solid castor and 24" PU rear wheel with aluminum frame



Model / **ALK868LB**

Name / Aluminum Manual Wheelchair

- 1) Nano Dryplating aluminum lightweight frame
- 2) Fixed armrest and footrest
- 3) 8"solid castor and 24" PU mag rear wheel



Model	Size (cm)																																																																																																									
-------	--------------	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	--	--	--	--	--	--	--	--	--	--



Model / **ALK954LBGC**
Name / Functional Aluminum Manual Wheelchair

- 1) Nano Dryplating aluminum lightweight frame
- 2) Adjustable reclining high back
- 3) Flip-up armrest, elevating long legrest
- 4) 8" solid front castor, 24" PU rear wheel with plastic handrim

Model	Size (cm)															
ALK954LBGC	62/67	35	41/46	59	20	49	93	118	42	46	100	86x35x88	18	20		



Model / **ALK908LAJP**
Name / Functional Aluminum Manual Wheelchair

- 1) Nano Dryplating aluminum lightweight frame
- 2) Drop back handle with united brake
- 3) Flip-up armrest, detachable footrest
- 4) 6" PVC front castor, 22" PU rear wheel with aluminum frame



Model / **ALK972LP-24"**
Name / Functional Aluminum Manual Wheelchair

- 1) Nano Dryplating, aluminum lightweight frame
- 2) Fixed armrest, detachable short footrest
- 3) 8" PVC front castor, 24" PU rear wheel with aluminum frame



Model / **ALK972LABJ-24"**
Name / Functional Aluminum Manual Wheelchair

- 1) Aluminum lightweight frame, double seat and back cushion (inner net seat for summer use and outer detachable soft cushion for winter use)
- 2) Fixed armrest and detachable footrest.
- 3) 8" PVC front castor and 24" PU mag rear wheel with plastic handrim
- 4) Drop back handle with united brake

Model	Size (cm)															
ALK908LAJP	58/63	27	41/46	56	15	50	89	107	40	39	100	78x27x74	13	15		
ALK972LP-24"	60/65	28	41/46	59	20	47	89	105	40	41	100	80x26x89	12	14		
ALK972LABJ-24"	62/67	30	41/46	60	20	47	89	105	40	41	100	80x33x72	12	14		



Name / Functional Aluminum Manual Wheelchair

- 1) Nano Drylating aluminum lightweight frame
- 2) Flip up and adjustable armrest
- 3) Detachable footrest
- 4) 8" PVC front castor, 24" quick release PU mag rear wheel



Name / Functional Aluminum Manual Wheelchair

- 1) Aluminum lightweight frame
- 2) elevating armrest and elevating long footrest
- 3) 8" PVC front castor, 24" quick release PU mag rear wheel
- 4) Double soft seat and back cushion

Name / Functional Aluminum Manual Wheelchair

- 1) Nano Dryplating, aluminum lightweight frame
- 2) Flip-up armrest, detachable short footrest
- 3) 8" PVC front castor, 24" quick release PU mag rear wheel

Name / Functional Aluminum Manual Wheelchair

- 1) Aluminum lightweight frame
- 2) Flip-up armrest, detachable short footrest
- 3) 8" PVC front castor, 24" quick release PU rear wheel

Name / Functional Aluminum Manual Wheelchair

- 1) Aluminum lightweight frame
- 2) Flip-up armrest, detachable footrest
- 3) 8" solid front castor, 24" PU mag rear wheel
- 4) Double soft seat and back cushion

Name / Functional Aluminum Manual Wheelchair

- 1) Aluminum lightweight frame
- 2) Flip-up armrest, elevating long footrest
- 3) 8" solid front castor, 24" PU mag rear wheel
- 4) Double soft seat and back cushion



Model	Size (cm)																																																																																																															
-------	--------------	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	--



Model \ Size (cm)													
	(cm)	(cm)	(cm)	(cm)	(cm)	(cm)	(cm)	(cm)	(cm)	(cm)	(cm)	(cm)	(cm)
ALK903LBQC	59/64	28	41/46	60	20	49	92	108	40	41	100	80x28x91	15 17
ALK953LB	60/65	28	41/46	60	20	50	89	109	40	40	100	80x28x91	16 18
ALK953LBC	60/65	28	41/46	60	20	50	89	109	40	40	100	80x28x91	16 18

Model / ALK809

Name / Steel Manual Wheelchair

- 1) Nano Dryplating steel frame
- 2) Fixed armrest and footrest
- 3) 8" PVC front castor,24" solid rear wheel



Model / ALK809

Name / Steel Manual Wheelchair

- 1) Chromed steel frame
- 2) Fixed armrest and footrest
- 3) 8" PVC front castor,24" solid rear wheel



Model / ALK809F

Name / Steel Manual Wheelchair

- 1) Chromed or powder coated steel frame
- 2) Fixed armrest and detachable footrest
- 3) 8"PVC front castor,24" solid rear wheel



Model	Size (cm)														
		60/65	23	41/46	59	20	49	87	104	40	40	100	91x21x87	15	17
ALK809		60/65	23	41/46	59	20	49	87	104	40	40	100	91x21x87	15	17
ALK809F		60/65	23	41/46	59	20	49	91	104	40	43	100	81x23x80	16	18



Model / ALK815

Name / Steel Manual Wheelchair

- 1) Chromed or powder coated steel frame
- 2) Fixed armrest and fixed footrest
- 3) 6"PVC front castor,24" solid rear wheel



Model / ALK815F

Name / Steel Manual Wheelchair

- 1) Chromed or powder coated steel frame
- 2) Fixed armrest and fixed footplate
- 3) 6"PVC front castor,24" solid rear wheel



Model	Size (cm)														
															
ALK815	63	23	43	54	15	46	84	98	40	41	100	89x23x85	13.5	15.5	
ALK815F	63	23	43	54	15	46	84	98	40	43	100	75x23x73	14	16	



Model / **ALK809B**

Name / Steel Manual Wheelchair

- 1) Nano Dryplating steel frame
- 2) Fixed armrest and footrest
- 3) 8" PVC front castor, 24" PU mag rear wheel



Model / **ALK809BJ**

Name / Steel Manual Wheelchair

- 1) Nano Dryplating steel frame with united brake
- 2) Fixed armrest and footrest
- 3) 8" PVC front castor, 24" PU mag rear wheel



Model / **ALK809Y**

Name / Steel Manual Wheelchair

- 1) Chromed steel frame with folding hard seat
- 2) Fixed armrest and footrest
- 3) 8" PVC front castor, 24" solid rear wheel



Model / **ALK874**

Name / Steel Manual Wheelchair

- 1) Nano Dryplating or chromed steel frame
- 2) Fixed armrest and footrest
- 3) 8" PVC front castor, 24" solid rear wheel



Model / **ALK874B**

Name / Steel Manual Wheelchair

- 1) Nano Dryplating or chromed steel frame
- 2) Fixed armrest and footrest
- 3) 8" PVC front castor, 24" PU mag rear wheel



Model	Size (cm)	Wheelchair width	Wheelchair depth	Wheelchair height	Wheelchair weight	Wheelchair length	Wheelchair width	Wheelchair height	Wheelchair weight	Wheelchair length	Wheelchair width	Wheelchair height	Wheelchair weight
ALK809B	62/67	24	41/46	60	20	49	87	107	40	40	100	94x23x88	15 17
ALK809BJ	62/67	24	41/46	60	20	49	87	107	40	40	100	94x23x88	15 17
ALK809Y	60/65	28	41/46	59	20	50	87	106	42	37	100	93x28x94	17 19

Model	Size (cm)	Wheelchair width	Wheelchair depth	Wheelchair height	Wheelchair weight	Wheelchair length	Wheelchair width	Wheelchair height	Wheelchair weight	Wheelchair length	Wheelchair width	Wheelchair height	Wheelchair weight
ALK874	61/66	26	41/46	59	20	48	90	105	40	42	110	95x25x90	16 18
ALK874B	63/68	28	41/46	60	20	48	90	105	40	42	110	95x27x90	16 18



Model / **ALK875**

Name / Steel Manual Wheelchair

- 1) Chromed steel frame
- 2) Fixed armrest and footrest
- 3) 8" PVC front castor, 24" solid rear wheel



Model / **ALK972**

Name / Functional Steel Manual Wheelchair

- 1) Chromed steel frame
- 2) Fixed armrest and detachable footrest
- 3) 8" PVC front castor, 24" solid rear wheel



Model / **ALK986AP**

Name / Functional Steel Manual Wheelchair

- 1) Powder coated steel frame with drop back handle
- 2) Fixed armrest, detachable footrest
- 3) 6" PVC front castor, 24" PU rear wheel



Model / **ALK904**

Name / Functional Steel Manual Wheelchair

- 1) Nano Dryplating or chromed steel frame
- 2) Flip-up armrest and detachable footrest
- 3) 8" PVC front castor, 24" solid rear wheel



Model / **ALK904BC**

Name / Functional Steel Manual Wheelchair

- 1) Nano Dryplating or chromed steel frame
- 2) Flip-up armrest and detachable elevating footrest
- 3) 8" PVC front castor, 24" PU mag rear wheel



Model	Size (cm)														
ALK875	60/65	23	41/46	59	20	49	87	106	40	40	100	93x21x87	15	17	
ALK972	60/65	23	41/46	59	20	49	87	102	40	41	100	80x23x87	17	19	
ALK986AP	57/62	26	41/46	59	15	50	94	103	40	45	100	77x26x75	14	16	

Model	Size (cm)												Hx GW	
														
ALK904	60/65	26	41/46	59	20	50	91	102	40	41	100	78x26x91	16 18	
ALK904BC	62/67	29	41/46	60	20	50	91	102	40	41	100	78x29x91	16 18	

Model / **ALK903**

Name / Functional Steel Manual Wheelchair

- 1) Powder coated steel frame
- 2) Detachable armrest and footrest
- 3) 8" PVC front castor, 24" solid rear wheel

Model / **ALK903B**

Name / Functional Steel Manual Wheelchair

- 1) Nano Dryplating or chromed steel frame
- 2) Detachable armrest and footrest
- 3) 8" PVC front castor, 24" PU mag rear wheel



Model / ALK903BC

Name / Functional Steel Manual Wheelchair

- 1) Nano Dryplating or chromed steel frame
- 2) Detachment armrest and elevating footrest
- 3) 8" PVC front castor, 24" PU mag rear wheel

Model \ Size (cm)													
ALK903	60/65	26	41/46	59	20	50	91	105	40	41	100	78x26x91	16 18
ALK903B	60/67	29	41/46	60	20	50	91	105	40	41	100	78x29x91	16 18
ALK903BC	62/67	29	41/46	60	20	50	91	105	40	41	100	78x29x91	16 18

Model / **ALK961 B**

Name / Functional Steel Manual Wheelchair

- 1) Powder coated steel frame
- 2) Flip up armrest, detachable footplate
- 3) 7"PVC front castor, 24"PU mag rear wheel



Model / ALK962B

Name / Functional Steel Manual Wheelchair

- 1) Powder coated steel frame
- 2) Flip up armrest, detachable footrest
- 3) 7" PVC front castor, 24" PU mag rear wheel

Model	Size (cm)													
ALK961B		63	28	46	59	17.5	52	93	100	40	41	120	87x28x93	17.5 19.5
ALK962B		63	28	46	59	17.5	52	93	105	40	41	120	76x28x93	17 19



Model / ALK905

Name / Functional Steel Manual Wheelchair

- 1) Nano Dryplating or chromed steel frame
- 2) Detachable armrest and footrest
- 3) 200*50 PVC front castor,24" solid rear wheel



Model / ALK905B

Name / Functional Steel Manual Wheelchair

- 1) Nano Dryplating or chromed steel frame
- 2) Detachable armrest and footrest
- 3) 200*50 PVC front castor,24" PU mag rear wheel



Model / ALK905C

Name / Functional Steel Manual Wheelchair

- 1) Nano Dryplating or chromed steel frame
- 2) Detachable armrest and elevating footrest
- 3) 200*50 PVC front castor,24" solid rear wheel

Model	Size (cm)														
ALK905	60/65	26	41/46	59	20	50	91	102	40	41	100	78x26x91	16	18	
ALK905B	62/67	29	41/46	60	20	50	91	102	40	41	100	78x29x91	16	18	
ALK905C	60/65	26	41/46	59	20	50	91	102	40	41	100	78x26x91	16	18	



Model / ALK951

Name / Functional Steel Manual Wheelchair

- 1) Powder coated steel frame
- 2) Flip-up armrest and detachable footrest
- 3) 200*50 PVC front castor,24" pneumatic rear wheel



Model / ALK951B

Name / Functional Steel Manual Wheelchair

- 1) Nano Dryplating or chromed steel frame
- 2) Flip-up armrest and detachable footrest
- 3) 200*50 PVC front castor,24" PU mag rear wheel



Model / ALK951BC

Name / Functional Steel Manual Wheelchair

- 1) Nano Dryplating or chromed steel frame
- 2) Flip-up armrest and detachable elevating footrest
- 3) 200*50 PVC front castor,24" PU mag rear wheel

Model	Size (cm)												H/W	GW
														
ALK951	71	28	56	61	20	52	93	107	45	42	125	81x28x93	18	20
ALK951B	73	30	56	60	20	52	93	107	45	42	125	81x30x93	18	20
ALK951BC	73	30	56	60	20	52	93	107	45	42	125	86x30x93	18	20



Model / **ALK601GC**

Name / Functional Commode Wheelchair

- 1) Nano Dryplating or chromed steel frame
- 2) Detachable armrest and elevating footrest
- 3) 8" PVC front castor, 24" solid rear wheel
- 4) Reclining high back

Model	Size (cm)													
ALK601GC	65	31	46	59	20	53	127	118	43	45	100	88x31x97	25	27



Model / ALK690GC

Name / Functional Commode Wheelchair

- 1) Powder coated steel frame
- 2) Fixed armrest and elevating footrest
- 3) 200*50"PVC fat wheel,24"pneumatic rear wheel

Model \ Size (cm)													
ALK690GC	67	30	44	60	20	52	122	124	45	45	100	89x30x94	19 21



Model / ALK603GC

Name / Functional Commode Wheelchair

- 1) Chromed steel frame
- 2) Detachable armrest and elevating footrest
- 3) 8" PVC front castor, 24" solid rear wheel
- 4) Reclining high back



Model / ALK605GC

Name / Functional Commode Wheelchair

- 1) Chromed steel frame
- 2) Flip-up armrest and elevating footrest
- 3) 8" PVC front castor, 24" solid rear wheel
- 4) Reclining high back



Model / ALK609GC

Name / Functional Commode Wheelchair

- 1) Chromed steel frame
- 2) Detachable armrest and elevating footrest
- 3) 8" PVC front castor, 24" solid rear wheel
- 4) Reclining high back



- 1) Nano Drylating aluminum lightweight frame
- 2) Detachable armrest and elevating footrest
- 3) 8" PVC front castor, 24" PU mag rear wheel
- 4) Reclining high back

Model / ALK608LBGC

Name / Functional Commode Wheelchair

Model	Size (mm)														
															
ALK603GC	65	31	46	59	20	53	127	118	43	45	100	88x31x95	25	27	
ALK605GC	65	31	46	59	20	53	127	118	43	44	100	90x31x97	25	27	
ALK609GC	65	31	46	59	20	56	127	118	43	44	100	90x31x97	24	26	

Model	Size (cm)													
ALK608LBGC	69	32	46	59	20	54	124	113	43	45	100	89×34×98	21	23

Model / **ALK681**

Name / Functional Commode Wheelchair

- 1) Nano Dryplating steel frame
- 2) Flip-up armrest and detachable footrest
- 3) 8" PVC front castor and 24" solid rear wheel



Model / ALK681B

Name / Functional Commode Wheelchair

- 1) Nano Drylating steel frame
- 2) Flip-up armrest and detachable footrest
- 3) 8" PVC front castor and 24" PU mag rear wheel



Model / ALK681BC

Name / Functional Commode Wheelchair

- 1) Nano Dryplating steel frame
- 2) Flip-up armrest and detachable elevating footrest
- 3) 8" PVC front castor and 24" PU mag rear wheel



Model / ALK609

Name / Functional Commode Wheelchair

- 1) Nano Drylating or chromed steel frame
- 2) Fixed armrest and footrest
- 3) 8" PVC front castor and 24" solid rear wheel



Model	Size (cm)													(kg)	(cm)	NW	GV
ALK681		65	26	46	59	20	54	91	105	43	36	100	78x32x91	20	22		
ALK681B		67	29	46	60	20	54	91	105	43	36	100	78x34x91	20	22		
ALK681BC		67	29	46	60	20	54	91	105	43	36	100	78x34x91	20	22		
ALK609		63	23	44	59	20	54	87	106	43	33	100	93x26x87	19	21		



Model / ALK608

Name / Commode Wheelchair

- 1) Nano Dryplating steel frame
- 2) Fixed armrest and footrest
- 3) 8" PVC front castor and 24" solid rear wheel

Model \ Size (cm)	63	23	44	59	20	52	87	106	43	36	100	93x26x87	18	20
ALK608														

Model / **ALK691-4"**
Name / Commode Wheelchair

- 1) Powder coated steel frame
- 2) Detachable armrest and footrest
- 3) 4" solid front castor, 4" rear wheel with brake



Model / **ALK691-5"**
Name / Commode Wheelchair

- 1) Powder coated steel frame
- 2) Detachable armrest and footrest
- 3) 5" solid front castor, 5" rear wheel with brake



Model / **ALK671S**
Name / Commode Wheelchair



Model / **ALK672S**
Name / Commode Wheelchair



Model / **ALK692**
Name / Commode Wheelchair

- 1) Chromed steel frame
- 2) Downstroke armrest and detachable footrest
- 3) 5" solid front castor, 5" rear wheel with brake



Model / **ALK693**
Name / Commode Wheelchair

- 1) Chromed steel frame
- 2) Downstroke armrest and detachable footrest
- 3) 5" solid front castor, 5" rear wheel with brake



Model / **ALK699L**
Name / Commode Wheelchair

- 1) Aluminum lightweight frame
- 2) Fixed armrest and flip-up footrest
- 3) 5" solid front castor, 5" rear wheel with brake



Model / **ALK698**
Name / Commode Wheelchair

- 1) Powder coated or chromed foldable steel frame
- 2) Downstroke armrest and detachable footrest
- 3) 5" solid front castor, 5" rear wheel with brake



Model	Size (cm)	Seat Width	Seat Depth	Seat Height	Wheel Width	Wheel Depth	Wheel Height	Wheel Diameter	Wheel Spacing	Wheel Base	Weight (kg)	Weight (cm)	NW	GW
ALK691-4"	53	-	43	10	10	50	91	93	43	37	100	55x53x49	13	15
ALK691-5"	53	-	43	12.5	12.5	52	95	93	43	37	100	55x53x52	13	15
ALK692	55	-	43	12.5	12.5	54	95	93	43	37	100	57x55x52	14	16
ALK693	55	-	43	12.5	12.5	54	99	93	43	42	100	68x44x57	26	28

Model	Size (cm)	Seat Width	Seat Depth	Seat Height	Wheel Width	Wheel Depth	Wheel Height	Wheel Diameter	Wheel Spacing	Wheel Base	Weight (kg)	Weight (cm)	NW	GW
ALK671S	34	6	34	36.5	-	-	36.5	39	37	40	-	71*46*36	12	
ALK672S	34.5	11.5	34.5	39.3	-	-	76.4	44	37	40	-	54*35*67	14	
ALK699L	54	-	38	12.5	12.5	54	89	55	44	35	100	80x38x55	14	16
ALK698	55	24	44	12.5	12.5	50	95	93	43	39	100	55x29x48	12	14

Model / **ALK613**
Name / Commode Wheelchair



Model / **ALK616L-3"**
Name / Commode Wheelchair



Model / **ALK616L**
Name / Commode Wheelchair



Model / **ALK611L**
Name / Commode Wheelchair



Model / **ALK612**
Name / Commode Wheelchair



Model / **ALK670S**
Name / Commode Wheelchair



Model / **ALK672**
Name / Commode Wheelchair



Model / **ALK690**
Name / Commode Wheelchair



Model / **ALK615**
Name / Commode Wheelchair



Model / **ALK675**
Name / Commode Wheelchair



Model	Size (mm)	Wheelchair		Walker		Rollator		Stairlift		Transfer		Weight	
	W	D	W	D	W	D	W	D	W	D	kg	lb	kg
ALK613	56	18	38	7.5	7.5	45-55	72-82	50	44	35	100	50x22x74	8 10
ALK616L-3"	56	18	38	7.5	7.5	45-55	78-88	50	44	35	100	50x22x74	5.5 7
ALK616L	56	18	38	-	-	45-55	78-88	50	44	35	100	50x22x74	5.5 7
ALK670S	61		50			42.5	90	56	43	51	125	78x37x63	5 11
ALK672	57		46			44	84	52	44	40	125	55x38.5x68	6.5 14
ALK690	58		47			52	92	82	44	40	125	54.5x19x76	89 9

Model \ Size (cm)								N.W (kg)	PCS CN	G.W (kg.)
ALK611L	46	55	45-55	78-88	44	35	58x46x63	4.5	2	11
ALK612	50	56	45-55	72-82	44	35	50x22x74	7	1	9
ALK615	48	44	48	82	43	34	50x25x82	9	1	11
ALK675	48	53	41	82	44	35	48x20x81	6	2	7.5



ALK410L



ALK401L



ALK402L



ALK403L



ALK406L



ALK408L



ALK405L



ALK405L-3"

Model	Size (cm)	Seat Width	Seat Depth	Seat Height	Backrest Height	Backrest Width	Backrest Depth	N.W (kg)	PCS CN	G.W (kg)
ALK410L	29	42	23	23	23	-	47x26x48	3	4	13
ALK401L	48	75	40-50	74-84	40	36	70x50x47	4	2	10
ALK402L	48	51	38-56	71-89	31	38	60x36x60	3	4	14
ALK403L	45	51	49	55	31	-	51x20x78	3	2	7
ALK406L	52	41	42-52	76-86	38	36	53x50x44	3	2	8
ALK408L	51	51	43-53	77-87	31	38	51x36x60	3	4	14
ALK405L	64	41	46	81	38	36	105x56x32	6	2	14
ALK405L-3"	64	41	56	91	38	36	109x56x32	7	2	16



ALK713L



ALK713LN



ALK715L



ALK762L



ALK762LN



ALK765L



ALK713L-5"



ALK765L-5"



ALK761L

Size (cm)						N.W (kg)	PCS CN	G.W (kg)
Model					(cm)			
ALK713L	58	10	47	78-95	54x30x80	2.3	4	12
ALK713LN	58	10	47	78-95	54x18x80	2.3	2	6
ALK715L	58	10	47	78-95	54x18x80	2.3	2	6
ALK762L	57	10	49	78-95	54x18x80	3.5	2	8
ALK762LN	55	10	50	74-92	54x18x77	3.5	2	8
ALK765L	59	12	45	78-95	54x18x80	2.5	2	6
ALK713L-5"	58	10	51	78-95	53x18x76	2.7	2	7
ALK765L-5"	59	12	52	80-97	53x18x76	3	2	7
ALK761L	58	10	48	75-92	54x19x78	3.25	2	7



ALK713L-S



ALK766L-S



ALK766L-S-3''



ALK766L-L



ALK763L-S



ALK763L-L

Model	Size (cm)					N.W (kg)	PCS CN	G.W (kg)
ALK713L-S	46	10	38	44-57	51x22x61	1.6	2	4
ALK766L-S	51	20	38	50-60	51x22x61	1.4	2	3
ALK766L-S-3''	51	20	38	50-60	51x19x61	1.4	2	3
ALK766L-L	51	13	53	83-93	75x22x95	1.8	2	4
ALK763L-S	57	16	73	47-65	53x30x58	4	2	10
ALK763L-L	65	18	84	77-97	67x34x91	5	2	11



ALK328



ALK327



ALK320L



ALK322L



ALK325L



ALK326L

Model	Size (cm)					N.W (kg)	PCS CN	G.W (kg)
ALK328	75	24	57	84-94	65x24x72	7	1	8.5
ALK327	62	24	68	88-98	68x24x60	7	1	9
ALK320L	61	22	60	80-91	61x22x83	6.5	1	8
ALK322L	61	20	64	75-86	61x20x80	6	1	7
ALK325L	61	18	72	82-93	54x20x54	6.5	1	8
ALK326L	61	22	70	80-91	61x22x86	7	1	8.5



ALK520LALK529LALK539L ALK530L ALK549L ALK528L ALK527L ALK536L



ALK525L (S) ALK525L (M) ALK525L (L)



ALK538LALK532LALK533L ALK537L ALK523L

Model	Size (cm)	Height(CM)		N.W (kg)	PCS CN	G.W (kg)
	(cm)					
ALK520L	71-94	62x15x28	0.35	20	8	
ALK529L	71-94	64x15x28	0.35	20	8	
ALK539L	68-91	62x13x28	0.35	20	8	
ALK530L	78-101	72x14x35	0.35	20	8	
ALK549L	68-91	60x13x28	0.35	20	8	
ALK528L	77-100	72x17x36	0.3	20	8	
ALK527L	82-92	61x16x28	0.33	20	8	
ALK536L	125	35x12x26	0.15	20	4	
ALK538L	80-103	80x19x32	0.33	20	9	
ALK532L	87-110	84x29x33	0.53	20	12	
ALK533L	85-107	84x29x29	0.53	20	12	
ALK537L	89-112	93x28x32	0.44	20	11	
ALK523L	86-109	86x26x29	0.5	20	11	
ALK520L-A1	71-94	62x15x28	0.35	20	8	
ALK527L-A1	82-92	61x16x28	0.35	20	8	
ALK540L	86	86x26x32	0.6	10	7	
ALK511L	86	85x30x40	0.82	10	10	
ALK525L(S)	95-116	101x33x29	0.74	20	16	
ALK525L(M)	114-135	119x33x29	0.82	20	17	
ALK525L(L)	135-155	140x33x29	0.9	20	19	
ALK534	68-91	77x34x34	0.77	10	9	
ALK526	73-96	75x33x34	0.74	10	9	
ALK524	63-86	74x30x34	0.8	10	9	
ALK521	63-86	97x34x44	0.86	10	10	
ALK522L	74-97	79x29x48	0.5	10	6	
ALK526L	73-96	75x33x34	0.5	10	6	



ALK520L-A1ALK527L-A1 ALK540L ALK511L



ALK534 ALK526 ALK524 ALK521 ALK522L ALK526L




ALK001
Urinal



ALK002
Urinal




ALK003
Urinal



(A1) 3" wheel
(A2) 3" lucency sport wheel
(A3) 5" wheel(M8 screw)
(A4) 5" wheel(608 Bearings)
(A5) 6" wheel
(A6) 8" wheel
(A7) 8" wheel(flat style)
(A8) 8" PVC fat wheel
(A9) 8" PU fat wheel
(B1) 6" castor with Al fork
(B2) 8" castor with Al fork

(B3) 8" castor with Steel fork
(B4) 5" grey plastic castor
(B5) 5" grey plastic castor with brake
(B6) 514screw medical PPR castor
(B7) 514screw medical PPR castor with brake
(B8) 5" wholly plastic castor
(B9) 5" wholly plastic castor with brake
(B10) Training castor
(C1) 12" PU wheel
(C2) 12" foaming wheel
(C3) 16" PU wheel


(C4) 20" PU wheel
(C5) 22" PU wheel
(C6) 24" PU wheel
(C7) 22" solid wheel
(C8) 24" solid wheel
(C9) 22" PU mag wheel(black tire)
(C10) 24" PU mag wheel(grey tire)
(C11) 24" solid wheel with plastic handrim
(C12) 24" PU wheel with drum brake
(C13) 20" PU wheel use in electric wheelchair
(C14) 24" PU wheel use in electric wheelchair



(D1) Walker foam grip
(D2) Crutch foam grip
(D3) Armrest pad
(D4) Crutch foam grip (L300)
(D5) Rollator back grip
(D6) Underarm crutch grip
(D7) Underarm shoulder pad
(D8) Walker stripe grip
(D9) Black leather short armrest pad
Black leather long armrest pad24" PU
Oxford cloth short armrest pad
Oxford cloth long armrest pad
PU short bent armrest pad
Plastic long armrest pad
Blow molding straight armrest pad
PU lengthen armrest pad
Ø 13 foot rubber
Ø 16 foot rubberwheelchair

(E3) Ø 19 foot rubber
(E4) Ø 19 special shape foot rubber
(E5) Ø 22 foot rubber
(E6) Ø 22 special shape foot rubber
(E7) Ø 25 foot rubber
(E8) Ø 28 foot rubber
(E9) Ø 28 big foot rubber
(E10) Ø 28 inclined foot rubber
(E11) Ø 19 circle big foot rubber
(E12) Ø 19 square big foot rubber
(F1) Flat crutch tips
(F2) 939 crutch tips
(F3) 929 crutch tips
(F4) Bow handle
(F5) 923 handle
(F6) 937 hand support
(F7) Disc grip
(F8) Ø 25 lead pipe

(G1) Ø 22 steel long footrest assembly
(G2) Ø 22 aluminum long footrest assembly
(G3) Ø 22 steel short footrest assembly
(G4) Ø 22 aluminum short footrest assembly
(G5) Ø 25 aluminum short footrest assembly
(G6) Aluminum pedals assembly
(G7) Plastic pedals assembly
(G8) Steel footrest assembly
(H1) 8" aluminum pedal
(H2) 7"aluminum pedal
(H3) Plastic pedal
(H4) 5" plastic pedal
(H5) Eccentricity pedal
(H6) 7"steel footrest
(H7) 6201 bearing
(H8) Bowl steel bearing
(H9) 6001 bearing
608 double convex bearing




- | | | | | | |
|-------|--|-------|---|------|---|
| (I1) | Steel brake assembly | (I13) | Ø 22 steel back loopers | (J7) | Aluminum flat armrest |
| (I2) | Steel united brake assembly | (I14) | Ø 22 aluminum back loopers | (J8) | Detachable aluminum short armrest |
| (I3) | A aluminum brake assembly | (I15) | Ø 25 aluminum back loopers | (J9) | Aluminum lift armrest assembly |
| (I4) | B aluminum brake assembly | (I16) | Drum brake assembly | (K1) | Toilet seatrest |
| (I5) | C aluminum brake assembly | (I17) | Steel drum assembly | (K2) | Blow seatrest |
| (I6) | D aluminum brake assembly | (I18) | Quick-release axle | (K3) | Commode seatrest |
| (I7) | Aluminum levers | (J1) | Detachable steel short armrest assembly | (K4) | High temperature pressure-shaped commode seatrest |
| (I8) | 8" steel levers | (J2) | Detachable steel armrest assembly | (K5) | Rollator seatrest |
| (I9) | 6" steel levers | (J3) | Detachable steel lengthen armrest assembly | (K6) | Round commode |
| (I10) | 8" aluminum brake for elec. wheelchair | (J4) | Detachable aluminum lengthen armrest assembly | (K7) | Bottomless commode |
| (I11) | 8" plastic levers | (J5) | Aluminum curved armrest | (K8) | Square commode |
| (I12) | 6" plastic levers | (J6) | Aluminum inclined armrest | (K9) | Square long commode |

Wheelchair upholstery



Nano Dryplating



Nano Dryplating has excellent corrosion resistance and adhesion. The brightness is more than ordinary powder and can achieve the effect of the chromed. Nano Dryplating can meet the latest EU environmental standards.